

Product datasheet

Anti-SHPRH antibody [2994C1a] ab84856

画像数 1

製品の概要

製品名	Anti-SHPRH antibody [2994C1a]
製品の詳細	Mouse monoclonal [2994C1a] to SHPRH
由来種	Mouse
アプリケーション	適用あり: WB, Dot blot
種交差性	交差種: Recombinant fragment 交差が予測される動物種: Human
免疫原	Recombinant fragment derived from the central region of Human SHPRH

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	2994C1a
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab84856** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution. Predicted molecular weight: 193 kDa.
Dot blot		Use at an assay dependent dilution.

ターゲット情報

機能	E3 ubiquitin-protein ligase involved in DNA repair. Upon genotoxic stress, accepts ubiquitin from the UBE2N-UBE2V2 E2 complex and transfers it to 'Lys-164' of PCNA which had been monoubiquitinated by UBE2A/B-RAD18, promoting the formation of non-canonical poly-ubiquitin chains linked through 'Lys-63'.
組織特異性	Broadly expressed.
パスウェイ	Protein modification; protein ubiquitination.
配列類似性	Belongs to the SNF2/RAD54 helicase family. Contains 1 H15 (linker histone H1/H5 globular) domain. Contains 1 helicase ATP-binding domain. Contains 1 helicase C-terminal domain. Contains 1 PHD-type zinc finger. Contains 1 RING-type zinc finger.
ドメイン	The RING finger mediates E3 ubiquitin ligase activity.

画像



Anti-SHPRH antibody [2994C1a] (ab84856) + immunogen
(recombinant fragment of Human SHPRH)

Predicted band size: 193 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors