

Product datasheet

Anti-SYTL3 antibody ab77020

画像数 1

製品の概要

製品名	Anti-SYTL3 antibody
製品の詳細	Mouse polyclonal to SYTL3
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	SYTL3 (AAI00831.1, 1 a.a. ~ 610 a.a) full-length human protein.
ポジティブ・コントロール	SYTL3 transfected 293T cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab77020** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 69 kDa.

ターゲット情報

機能	May act as Rab effector protein and play a role in vesicle trafficking. Binds phospholipids in the presence of calcium ions.
配列類似性	Contains 2 C2 domains. Contains 1 RabBD (Rab-binding) domain.
細胞内局在	Endomembrane system.

画像



All lanes : Anti-SYTL3 antibody (ab77020) at 5 µg/ml

Lane 1 : SYTL3 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size: 69 kDa

Observed band size: 67 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors