

Product datasheet

Anti-Salmonella A+G serotype antibody ab78686

製品の概要

製品名	Anti-Salmonella A+G serotype antibody
製品の詳細	Rabbit polyclonal to Salmonella A+G serotype
由来種	Rabbit
特異性	ab78686 recognizes somatic O antigens 2, 3, 4, 5, 7, 8, 9, 10, 12, 15, 19, 34, 46, 11 and 13.
アプリケーション	適用あり: Agglutination
種交差性	Reacts with Salmonella
免疫原	Tissue/ cell preparation of Salmonella (heat inactivated).
ポジティブ・コントロール	Antigenic suspension.
特記事項	Reagents should be allowed to stand at 15°C - 25°C for at least 30 minutes prior to use.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.08% Sodium Azide Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab78686** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Agglutination		Use at an assay dependent concentration. See attached protocol.

ターゲット情報

関連性

Salmonella is a genus of rod-shaped, Gram-negative, non-spore forming, predominantly motile enterobacteria which cause many types of infections from mild gastroenteritis to life threatening typhoid fever. *Salmonella* can be serotyped according to somatic antigens O (group O: O2-O67 groups), surface Vi, phase I flagellar and phase II flagellar (H types: H-a- H-z81) antigens. The somatic O antigen are heat stable and upon which grouping of the organism is based. The H antigen are heat labile and usually associated with motility. H antigen exist in 2 phases, phase I and phase II.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors