

Product datasheet

Anti-Sca1 / Ly6A/E antibody [D7] (PerCP/Cy5.5®) ab93538

画像数 1

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-Sca1 / Ly6A/E antibody [D7] (PerCP/Cy5.5®) |
| 製品の詳細 | Rat monoclonal [D7] to Sca1 / Ly6A/E (PerCP/Cy5.5®) |
| 由来種 | Rat |
| 標識 | PerCP/Cy5.5®. Ex: 482nm, Em: 690nm |
| アプリケーション | 適用あり: Flow Cyt |
| 種交差性 | 交差種: Mouse |
| 免疫原 | IL2 dependent T cell line CTL-L (Mouse). |
| ポジティブ・コントロール | C57Bl/6 Mouse splenocytes |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | pH: 7.20 Preservative: 0.09% Sodium azide Constituents: 0.87% Sodium chloride, 0.1% Gelatin, PBS |
| 精製度 | Protein G purified |
| ポリ/モノ | モノクローナル |
| クローン名 | D7 |
| アイソタイプ | IgG2a |
| 軽鎖の種類 | kappa |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab93538** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
| Flow Cyt | | |

追加情報

Flow Cyt: Use 0.125 - 0.25 µg for 10^5 - 10^8 cells in a 100 µl volume.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

| | |
|-------|--|
| 機能 | T-cell activation. |
| 組織特異性 | Widely expressed. |
| 配列類似性 | Contains 1 UPAR/Ly6 domain. |
| 翻訳後修飾 | O-glycosylated. Not N-glycosylated. Not phosphorylated. |
| 細胞内局在 | Cell membrane. |

画像



Flow cytometry staining of C57Bl/6 splenocytes with 0.25 µg of PerCP-Cy5.5 Rat IgG2a Isotype Control (open histogram) or 0.125 µg of ab93538 (colored histogram). Cells in the lymphocyte gate were used for analysis.

Flow Cytometry - Anti-Sca1 / Ly6A/E antibody [D7]
(PerCP/Cy5.5®) (ab93538)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors