


Product datasheet

Anti-SerpinB2 / PAI-2 antibody ab105739

画像数 1

製品の概要

製品名	Anti-SerpinB2 / PAI-2 antibody
製品の詳細	Rabbit polyclonal to SerpinB2 / PAI-2
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human 交差が予測される動物種: Monkey 
免疫原	Synthetic peptide corresponding to Human SerpinB2/ PAI-2 (internal sequence). Database link: NP_002566.1
ポジティブ・コントロール	HeLa whole cell lysates.
特記事項	Expiration date is one (1) year from date of opening product. Previously labelled as SerpinB2.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.88% Sodium chloride
精製度	Protein A purified
特記事項 (精製)	ab105739 is affinity purified from monospecific antiserum by immunoaffinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab105739** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/2000. Detects a band of approximately 42 kDa (predicted molecular weight: 47 kDa).
ELISA		1/100000 - 1/150000.

ターゲット情報

機能	Inhibits urokinase-type plasminogen activator. The monocyte derived PAI-2 is distinct from the endothelial cell-derived PAI-1.
配列類似性	Belongs to the serpin family. Ov-serpin subfamily.
翻訳後修飾	The signal sequence is not cleaved.
細胞内局在	Cytoplasm. Secreted > extracellular space.

画像



Anti-SerpinB2 / PAI-2 antibody (ab105739) at 1/500 dilution + HeLa whole cell lysates at 50 µg

Secondary

Goat-anti-Rabbit HRP at 1/5000 dilution

Developed using the ECL technique.

Predicted band size: 47 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors