

## Product datasheet

### Anti-TFIIB antibody [EPR4589] ab109106

リコンビナント RabMAb®

1 References 画像数 1

#### 製品の概要

製品名	Anti-TFIIB antibody [EPR4589]
製品の詳細	Rabbit monoclonal [EPR4589] to TFIIB
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB, IP <b>適用なし:</b> Flow Cyt, ICC or IHC-P
種交差性	<b>交差種:</b> Mouse, Rat, Human
免疫原	Synthetic peptide within Human TFIIB aa 300-400 (C terminal). The exact sequence is proprietary.
ポジティブ・コントロール	HepG2, HeLa, Jurkat, K562, C6, RAW 264.7, and PC12 cell lysates.
特記事項	

Our RabMAb® technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb® patents](#)

This product is a recombinant rabbit monoclonal antibody.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	pH: 7.20 Preservative: 0.05% Sodium azide Constituents: 0.1% BSA, 40% Glycerol, 9.85% Tris glycine, 50% Tissue culture supernatant
精製度	Tissue culture supernatant
ポリ/モノ	モノクローナル
クローン名	EPR4589
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab109106** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/10000. Predicted molecular weight: 35 kDa.
IP		1/10 - 1/100.

**追加情報** Is unsuitable for Flow Cyt, ICC or IHC-P.


### ターゲット情報

**機能** General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.

**配列類似性** Belongs to the TFIIIB family.  
Contains 1 TFIIIB-type zinc finger.

**細胞内局在** Nucleus.

### 画像



**All lanes :** Anti-TFIIIB antibody [EPR4589] (ab109106) at 1/1000 dilution

- Lane 1 :** HepG2 cell lysate
- Lane 2 :** HeLa cell lysate
- Lane 3 :** Jurkat cell lysate
- Lane 4 :** K562 cell lysate
- Lane 5 :** C6 cell lysate
- Lane 6 :** RAW 264.7 cell lysate
- Lane 7 :** PC12 cell lysate

Lysates/proteins at 10 µg per lane.

**Predicted band size:** 35 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors