

## Product datasheet

# Anti-TMEM143 antibody ab121431

画像数 3

### 製品の概要

製品名	Anti-TMEM143 antibody
製品の詳細	Rabbit polyclonal to TMEM143
由来種	Rabbit
アプリケーション	<b>適用あり:</b> ICC/IF, IHC-P, WB
種交差性	<b>交差種:</b> Human
免疫原	fragment, corresponding to internal sequence amino acids 57-197 of Human TMEM143, isoform 1.
ポジティブ・コントロール	Human colon tissue; Transfected HEK293T cells.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab121431** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		1/250 - 1/500.

## ターゲット情報


細胞内局在 Membrane.

## 画像




Immunofluorescent staining of human cell line U-251 MG shows positivity in nucleoli.

Immunocytochemistry/ Immunofluorescence - Anti-TMEM143 antibody (ab121431)



ab121431 at a 1/250 dilution staining  
TMEM43 in paraffin-embedded Human colon  
tissue, by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-  
embedded sections) - Anti-TMEM143 antibody  
(ab121431)



**All lanes :** Anti-TMEM143 antibody  
(ab121431) at 1/250 dilution

**Lane 1 :** Vector only transfected HEK293T  
lysate

**Lane 2 :** TMEM143-transfected HEK293T  
over-expression Lysate (Co-expressed with a  
C-terminal myc-DDK tag ~3.1 kDa)

Western blot - Anti-TMEM143 antibody (ab121431)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors