


Product datasheet

Anti-TMEM147 antibody ab97624

画像数 4

製品の概要

製品名	Anti-TMEM147 antibody
製品の詳細	Rabbit polyclonal to TMEM147
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P, ICC/IF
種交差性	交差種: Mouse, Human 交差が予測される動物種: Rat, Cow, Zebrafish 
免疫原	Synthetic peptide corresponding to a region within amino acids 87 and 101 of Human TMEM147.
ポジティブ・コントロール	WB: Raji or NIH 3T3 cell lysate; ICC/IF: A431 cells; IHC-P: NCI-N87 xenograft.

法規制情報

医薬用外毒物

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab97624** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 25 kDa.
IHC-P		1/100 - 1/500.
ICC/IF		1/100 - 1/200.

## ターゲット情報

細胞内局在 Endoplasmic reticulum membrane.


## 画像



Anti-TMEM147 antibody (ab97624) at 1/1500 dilution + Raji whole cell lysate at 30  $\mu$ g

**Predicted band size:** 25 kDa


12% SDS PAGE



Anti-TMEM147 antibody (ab97624) at 1/1000 dilution + NIH 3T3 whole cell lysate at 30  $\mu$ g

**Predicted band size:** 25 kDa

12% SDS PAGE




1/200 ab97624 staining TMEM147 in paraformaldehyde-fixed A431 cells by immunofluorescence.

The lower image was co-stained with Hoechst 33342.



Immunocytochemistry/ Immunofluorescence - Anti-TMEM147 antibody (ab97624)



1/500 ab97624 staining TMEM147 in paraffin-embedded NCI-N87 xenograft by immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TMEM147 antibody (ab97624)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors