

Product datasheet

Anti-UBQLN2 antibody ab108059

画像数 1

製品の概要

製品名	Anti-UBQLN2 antibody
製品の詳細	Rabbit polyclonal to UBQLN2
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Guinea pig, Cow, Cat, Dog 
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 81-130 (LKDQDTLIQH GIHDGLTVHL VIKSQNRPQG QSTQPSNAAG TNTTSASTPR) of Human UBQLN2 (NP_038472). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	HepG2 cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab108059** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 66 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

関連性	UBQLN2 (Ubiquilin 2) is a ubiquitin-like protein which contains an N-terminal ubiquitin-like (UBL) domain and a C-terminal ubiquitin-associated (UBA) domain. Ubiquilins interact with both proteasomes and ubiquitin ligases, and are thought to link the ubiquitination machinery to the proteasome to affect in vivo protein degradation. UBQLN2 increases the half-life of proteins destined to be degraded by the proteasome.
細胞内局在	Cytoplasmic and Nuclear. Note=Where it colocalizes with the proteasome. Associated with fibers in mitotic cells.

画像



Anti-UBQLN2 antibody (ab108059) at 1 µg/ml
+ HepG2 cell lysate at 10 µg

Predicted band size: 66 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors