

## Product datasheet

# Anti-UTP3 antibody [2950C2a] ab60063

### 画像数 1

#### 製品の概要

---

|          |  |
|----------|--|
| 製品名      | Anti-UTP3 antibody [2950C2a]   |
| 製品の詳細    | Mouse monoclonal [2950C2a] to UTP3   |
| 由来種      | Mouse  |
| アプリケーション | <b>適用あり:</b> WB, Dot blot  |
| 種交差性     | <b>交差種:</b> Human  |
| 免疫原      | Recombinant fragment, sequence taken from near the C terminus of Human UTP3. |

#### 製品の特性

---

|          |  |
|----------|--|
| 製品の状態    | Liquid   |
| 保存方法     | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.   |
| バッファー    | pH: 7.40<br>Preservative: 0.05% Sodium azide<br>Constituents: 0.0204% Monobasic dihydrogen potassium phosphate, 0.1136% Dibasic monohydrogen sodium phosphate, 0.0225% Potassium chloride, 0.03% Potassium phosphate, 0.1312% Sodium phosphate, 0.812% Sodium chloride, 1% BSA |
| 精製度      | Protein G purified   |
| 特記事項(精製) | ab60063 was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum and filtered through a 0.22µm membrane.  |
| ポリモノ     | モノクローナル  |
| クローン名    | 2950C2a  |
| アイソタイプ   | IgG1   |

#### アプリケーション

---

Our [Abpromise guarantee](#) covers the use of **ab60063** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB

Dot blot

## 追加情報

Dot: Use at an assay dependent dilution.

WB: Use at an assay dependent dilution. Detects a band of approximately 39 kDa (fragment), (predicted molecular weight: 55 kDa).

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

## 関連性

UTP3, also known as SAS10 (Something about silencing protein 10) is essential for gene silencing; having a role in the structure of silenced chromatin. It also plays a role in the developing brain.

## 細胞内局在

Nuclear

## 画像



Western blot - Anti-UTP3 antibody [2950C2a] (ab60063)

Anti-UTP3 antibody [2950C2a] (ab60063) + immunising recombinant protein

**Predicted band size:** 55 kDa

**Observed band size:** 39 kDa

The molecular weight of the band on the western blot does not correspond to the molecular weight of the natural protein because only a fragment of the protein was used.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors