

Product datasheet

Anti-WDR73 antibody ab103864

1 References 画像数 2

製品の概要

製品名	Anti-WDR73 antibody
製品の詳細	Rabbit polyclonal to WDR73
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 52-82 of Human WDR73.
ポジティブ・コントロール	MK562 cell lysate; Human brain tissue

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
特記事項(精製)	ab103864 is purified through a protein A column, followed by peptide affinity purification.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab103864** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------


WB	1/100 - 1/500. Predicted molecular weight: 42 kDa.
----	--

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報


関連性 WDR73 belongs to the WD repeat WDR73 family and contains 3 WD repeats.

画像



Anti-WDR73 antibody (ab103864) at 1/100 dilution + MK562 cell lysate at 35 µg

Predicted band size: 42 kDa



ab103864, at 1/50 dilution, staining WDR73 in formalin-fixed, paraffin-embedded Human brain tissue by Immunohistochemistry using a peroxidase conjugated secondary antibody and DAB staining.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors