


## Product datasheet

# Anti-WW domain binding protein 4 antibody ab108144

### 画像数 1

#### 製品の概要

製品名	Anti-WW domain binding protein 4 antibody
製品の詳細	Rabbit polyclonal to WW domain binding protein 4
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Zebrafish 
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MADYWKSQPK KFCDYCKCWI ADNRPVSVEFH ERGKNHKENV AKRISEIKQK) of Human WW domain binding protein 4 (NP_009118). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	721_B cell lysate

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab108144** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 43 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## ターゲット情報

機能	Promotes pre-mRNA splicing. A spliceosome-associated protein; may play a role in cross-intron bridging of U1 and U2 snRNPs in the mammalian A complex.
配列類似性	Contains 1 matrin-type zinc finger. Contains 2 WW domains.
ドメイン	The WW domain recognizes the proline, glycine and methionine-rich (PGM) motif present in the splicing factors, as well as the Arg/Gly-rich-flanked Pro-rich domains found in several WW domain-binding proteins.
細胞内局在	Nucleus speckle.

## 画像



Anti-WW domain binding protein 4 antibody (ab108144) at 1 µg/ml  
+ 721\_B cell lysate at 10 µg

**Predicted band size:** 43 kDa

Gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors