

## Product datasheet

# Anti-ZBTB48 antibody ab84023

### 画像数 1

#### 製品の概要

製品名	Anti-ZBTB48 antibody
製品の詳細	Rabbit polyclonal to ZBTB48
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide within residues: EPAGLEEEEV SRLTGLVPRD QEPRGSHSPQ RPQLHSPAQS EGPSSLCGKL , corresponding to N terminal amino acids 144-193 of Human ZBTB48 (NP_005332) <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	293T cell lysate.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab84023** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項


アプリ ケーショ ン	Abreviews	特記事項
------------------	-----------	------

WB Use a concentration of 1 µg/ml. Detects a band of approximately 77 kDa (predicted molecular weight: 77 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## ターゲット情報

機能	Binds to and regulates the J and/or S elements in MHC II promoter.
組織特異性	Detected in adrenal gland and neuroblastoma.
配列類似性	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 1 BTB (POZ) domain. Contains 11 C2H2-type zinc fingers.
細胞内局在	Nucleus.

## 画像



Anti-ZBTB48 antibody (ab84023) at 1 µg/ml + 293T cell lysate at 10 µg

**Predicted band size:** 77 kDa

**Observed band size:** 77 kDa

**Additional bands at:** 30 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors