

Product datasheet

Anti-ZC3HAV1L antibody ab121624

画像数 2

製品の概要

製品名	Anti-ZC3HAV1L antibody
製品の詳細	Rabbit polyclonal to ZC3HAV1L
由来種	Rabbit
アプリケーション	適用あり: ICC/IF, IHC-P
種交差性	交差種: Human
免疫原	antigen sequence: CKRSHQLIHA ASLKLLQDQG LNIPSVVNFQ IISTYKHMKL HKMLENTDNS SPSTHSQGL EKQGVHAAGA AEAGPLASVP AQSAAKPCPV, corresponding to amino acids 207-296 of Human ZC3HAV1L. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	Human pancreas tissue.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab121624** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項


アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment
IHC-P		1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. Use buffer at pH 6 for antigen retrieval.

## ターゲット情報

### 翻訳後修飾


Phosphorylated upon DNA damage, probably by ATM or ATR.

## 画像



Immunofluorescent staining of Human cell line U-251MG shows positivity in cytoplasm. Recommended concentration of ab121624 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-ZC3HAV1L antibody (ab121624)



ab121624 at 1/500 staining ZC3HAV1L in Paraffin Embedded Human pancreas tissue by Immunohistochemistry.  
Note: strong cytoplasmic positivity in exocrine glandular cells.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ZC3HAV1L antibody (ab121624)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors