

Product datasheet

Anti-ZMYM2 antibody ab106624

画像数 3

製品の概要

製品名	Anti-ZMYM2 antibody
製品の詳細	Rabbit polyclonal to ZMYM2
由来種	Rabbit
特異性	At least three isoforms of ZMYM2 are known to exist; ab106624 will detect only the longest isoform. It is predicted to not cross react with other ZMYM protein family members.
アプリケーション	適用あり: ICC/IF, WB
種交差性	交差種: Mouse, Rat, Human
免疫原	A 17 amino acid peptide near the amino terminus of Human ZMYM2 (CAH71822).
ポジティブ・コントロール	EL4 cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (stable for up to 12 months).
バッファー	Preservative: 0.02% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab106624** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 10 µg/ml.
WB		Use a concentration of 0.125 - 0.25 µg/ml. Predicted molecular weight: 155 kDa.

ターゲット情報


機能	May function as a transcription factor.
関連疾患	Note=A chromosomal aberration involving ZMYM2 may be a cause of stem cell leukemia lymphoma syndrome (SCLL). Translocation t(8;13)(p11;q12) with FGFR1. SCLL usually presents as lymphoblastic lymphoma in association with a myeloproliferative disorder, often accompanied by pronounced peripheral eosinophilia and/or prominent eosinophilic infiltrates in the affected bone marrow.
配列類似性	Contains 9 MYM-type zinc fingers.
細胞内局在	Nucleus.
製品の状態	There are 2 isoforms produced by alternative splicing.

画像



Immunofluorescence of ZMYM2 in HeLa cells with ab106624 at 20 ug/mL.

Immunocytochemistry/ Immunofluorescence - Anti-ZMYM2 antibody (ab106624)



Lane 1 : Anti-ZMYM2 antibody (ab106624) at 0.125 µg/ml

Lane 2 : Anti-ZMYM2 antibody (ab106624) at 0.25 µg/ml

All lanes : EL4 cell lysate

Lysates/proteins at 15 µg per lane.

Predicted band size: 155 kDa

Western blot - Anti-ZMYM2 antibody (ab106624)



Immunocytochemistry of ZMYM2 in HeLa cells
with ab106624 at 10 ug/mL.

Immunocytochemistry - Anti-ZMYM2 antibody
(ab106624)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors