

## Product datasheet

# Anti-ZNF519 antibody - N-terminal ab91248

### 画像数 1

#### 製品の概要

|              |   |
|--------------|---|
| 製品名          | Anti-ZNF519 antibody - N-terminal   |
| 製品の詳細        | Rabbit polyclonal to ZNF519 - N-terminal  |
| 由来種          | Rabbit  |
| アプリケーション     | <b>適用あり:</b> ELISA, WB  |
| 種交差性         | <b>交差種:</b> Human   |
| 免疫原          | KLH conjugated synthetic peptide selected from the N terminal region of Human ZNF519. |
| ポジティブ・コントロール | T47D cell line lysates  |

#### 製品の特性

|           |   |
|-----------|---|
| 製品の状態     | Liquid  |
| 保存方法      | Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.        |
| バッファー     | Preservative: 0.09% Sodium Azide<br>Constituents: PBS                           |
| 精製度       | Immunogen affinity purified   |
| 特記事項 (精製) | Purified through a protein A column, followed by peptide affinity purification. |
| ポリモノ      | ポリクローナル   |
| アイソタイプ    | IgG   |

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab91248** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
| ELISA    |           |      |
| WB       |           |      |

追加情報 ELISA: 1/1000.

WB: 1/100 - 1/500. Predicted molecular weight: 63 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報


---

**関連性** ZNF519 may be involved in transcriptional regulation.

**細胞内局在** Nuclear

## 画像

---



Anti-ZNF519 antibody - N-terminal (ab91248)  
at 1/100 dilution + T47D cell line lysate at 35  
µg

**Predicted band size: 63 kDa**

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors