

Product datasheet

Anti-Parvin alpha antibody ab102646

画像数 2

製品の概要

製品名	Anti-Parvin alpha antibody
製品の詳細	Rabbit polyclonal to Parvin alpha
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Mouse, Human
免疫原	Synthetic peptide conjugated to KLH, between C-terminal amino acids 330-360 of Human Parvin alpha.
ポジティブ・コントロール	Human Skeletal Muscle tissue and Mouse Bladder Tissue lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項(精製)	ab102646 is purified through a protein A column, followed by peptide affinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab102646** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/100 - 1/500. Predicted molecular weight: 42 kDa.

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報

機能	Probably plays a role in the regulation of cell adhesion and cytoskeleton organization. Plays a role in ciliogenesis.
組織特異性	Widely expressed, with highest levels in heart, skeletal muscle, kidney and liver.
配列類似性	Belongs to the parvin family. Contains 2 CH (calponin-homology) domains.
細胞内局在	Cell junction > focal adhesion. Cell membrane. Cytoplasm > cytoskeleton. Constituent of focal adhesions.

画像



Western blot - Anti-Parvin alpha antibody (ab102646)

Anti-Parvin alpha antibody (ab102646) at
1/100 dilution + Mouse bladder tissue lysates
at 35 µg

Predicted band size: 42 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Parvin alpha antibody (ab102646)

IHC analysis in formalin fixed and paraffin embedded Human skeletal muscle tissue using [ab102644](#) at 1/100 followed by peroxidase conjugation of the secondary antibody and DAB staining.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors