

## Product datasheet

# Anti-Cotinine antibody [B948M] ab104810

### 製品の概要

製品名	Anti-Cotinine antibody [B948M]
製品の詳細	Mouse monoclonal [B948M] to Cotinine
由来種	Mouse
アプリケーション	<b>適用あり:</b> ELISA
免疫原	Cotinine-BSA

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 0.85% Sodium chloride, PBS, pH 7.2
精製度	DEAE-Chromatography
ポリ/モノ	モノクローナル
クローン名	B948M
アイソタイプ	IgG1

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab104810** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.

### ターゲット情報

**関連性** Cotinine is a chemical derived from the catabolism of nicotine, which is found in cigarette smoke. Since cotinine can be made only from nicotine, and since nicotine enters the body with

cigarette smoke, cotinine measurements can show how much cigarette smoke enters the body in a given time. Cotinine itself is not harmful.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors