

## Product datasheet

# Anti-Trichomonas foetus antibody [B1212M] ab105854

### 製品の概要

製品名	Anti-Trichomonas foetus antibody [B1212M]
製品の詳細	Mouse monoclonal [B1212M] to Trichomonas foetus
由来種	Mouse
特異性	Does not cross react with Trichomonas vaginalis.
アプリケーション	<b>適用あり:</b> ELISA, ICC/IF, Sandwich ELISA
種交差性	<b>交差種:</b> Other species
免疫原	Disrupted Trichomonas foetus

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.1% Sodium Azide Constituents: 0.01M PBS, pH 7.2
精製度	Protein A purified
特記事項(精製)	ab105854 is >90% pure
ポリ/モノ	モノクローナル
クローン名	B1212M
アイソタイプ	IgG2b

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab105854** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.

## アプリケーション

## Abreviews

## 特記事項

Sandwich ELISA

Use at an assay dependent dilution.

## ターゲット情報

### 関連性

*Tritrichomonas foetus* is a microscopic single-celled flagellated protozoan parasite that has traditionally been identified as a cause of reproductive disease in cattle (infertility, abortion and endometritis). It has been found all over the world, but the widespread use of artificial insemination in breeding cattle has led to the virtual elimination of this organism from the cattle population in many countries including the UK and much of Europe. *Tritrichomonas foetus* may also be an important cause of diarrhoea in cats. It can infect and colonises the large intestine, and can cause prolonged and intractable diarrhoea.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors