

Product datasheet

Anti-FKBP11 antibody [5] ab106198

画像数 1

製品の概要

製品名	Anti-FKBP11 antibody [5]
製品の詳細	Mouse monoclonal [5] to FKBP11
由来種	Mouse
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human
免疫原	Recombinant full length Human FKBP11 protein (NP_057678.1) expressed in E. coli system.
特記事項	Affinity = 7.68×10^9

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.02% Proclin Constituents: Tris HCl, Sodium citrate, pH 7.0
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	5
アイソタイプ	IgG2b

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab106198** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution. Predicted molecular weight: 22 kDa.
ELISA		Use at an assay dependent dilution.

ターゲット情報

機能	PPlases accelerate the folding of proteins during protein synthesis.
配列類似性	Belongs to the FKBP-type PPlase family. Contains 1 PPlase FKBP-type domain.
細胞内局在	Membrane.

画像



Lane 1 : Anti-FKBP11 antibody [5]
(ab106198) at 1/1000 dilution

Lane 2 : N/A

Lane 1 : Jurkat whole cell lysate.

Lane 2 : Molecular weight ladder.

Predicted band size: 22 kDa

Observed band size: 22 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors