


Product datasheet

Anti-TMBIM4 antibody ab110376

画像数 1

製品の概要

製品名	Anti-TMBIM4 antibody
製品の詳細	Rabbit polyclonal to TMBIM4
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Rat 
免疫原	Synthetic peptide derived from an internal sequence of Human TMBIM4.
ポジティブ・コントロール	HepG2 cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: PBS, 0.88% Sodium chloride, 50% Glycerol Note: without Mg ²⁺ and Ca ²⁺
精製度	Immunogen affinity purified
特記事項(精製)	ab110376 was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab110376** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション

Abreviews

特記事項


WB 1/500 - 1/1000. Predicted molecular weight: 27 kDa.

ターゲット情報

配列類似性 Belongs to the BI1 family.

細胞内局在 Membrane.

画像



All lanes : Anti-TMBIM4 antibody (ab110376)
at 1/500 dilution

Lane 1 : HepG2 cell lysate

Lane 2 : HepG2 cell lysate with immunizing
peptide at 10 µg

Lysates/proteins at 30 µg per lane.

Predicted band size: 27 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors