


Product datasheet

Anti-OR10H3 + OR10H4 antibody ab110436

画像数 1

製品の概要

製品名	Anti-OR10H3 + OR10H4 antibody
製品の詳細	Rabbit polyclonal to OR10H3 + OR10H4
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse 交差が予測される動物種: Human 
免疫原	Synthetic peptide derived from Human OR10H3 + OR10H4.
ポジティブ・コントロール	Extracts from NIH-3T3 cells

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: PBS, 0.88% Sodium chloride, 50% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab110436** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Predicted molecular weight: 36 kDa.

<b>関連性</b>	OR10H3 and OR10H4 are odorant receptors and belong to the G-protein coupled receptor 1 family.
<b>細胞内局在</b>	Cell membrane; Multi-pass membrane protein.

## 画像



**All lanes :** Anti-OR10H3 + OR10H4 antibody (ab110436) at 1/500 dilution

**Lane 1 :** Extracts from NIH-3T3 cells

**Lane 2 :** Extracts from NIH-3T3 cells with immunising peptide at 10 µg

Lysates/proteins at 30 µg per lane.

**Predicted band size:** 36 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors