

## Product datasheet

# Anti-TRIP13 antibody ab113008

### 画像数 1

#### 製品の概要

製品名	Anti-TRIP13 antibody
製品の詳細	Rabbit polyclonal to TRIP13
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 370-400 of Human TRIP13(NP_004228.1).
ポジティブ・コントロール	293 cell lysates nontransfected or transiently transfected with the TRIP13 gene.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium azide Constituent: 99% PBS
精製度	Ammonium Sulphate Precipitation
特記事項(精製)	ab113008 was prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab113008** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/50 - 1/100. Predicted molecular weight: 49 kDa.

## ターゲット情報

### 配列類似性

Belongs to the AAA ATPase family.

### 画像



**All lanes :** Anti-TRIP13 antibody (ab113008)  
at 1/50 dilution

**Lane 1 :** non-transfected 293 cell lysate

**Lane 2 :** TRIP13 transiently transfected 293  
cell lysate

Lysates/proteins at 2 µg per lane.

**Predicted band size:** 49 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors