

## Product datasheet

# Anti-CIP29 antibody ab115175

### 画像数 1

#### 製品の概要

製品名	Anti-CIP29 antibody
製品の詳細	Rabbit polyclonal to CIP29
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IP <b>適用なし:</b> WB
種交差性	<b>交差種:</b> Human <b>交差が予測される動物種:</b> Rat, Cow, Chimpanzee, Lizard, Rhesus monkey, Gorilla, Chinese hamster, Orangutan 
免疫原	Synthetic peptide corresponding to a region between residue 1 and 50 of Human CIP29 (NP_149073.1).
ポジティブ・コントロール	Whole cell lysate from HeLa cells.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.50 Preservative: 0.09% Sodium azide Constituent: 99% Tris citrate/phosphate
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab115175** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IP		Use at 2-10 µg/mg of lysate.

## 追加情報

Is unsuitable for WB.

## ターゲット情報

### 機能

Binds both single-stranded and double-stranded DNA with higher affinity for the single-stranded form. Specifically binds to scaffold/matrix attachment region DNA. Also binds single-stranded RNA. Enhances RNA unwinding activity of DDX39A. May participate in important transcriptional or translational control of cell growth, metabolism and carcinogenesis.

### 組織特異性

Low expression in spleen, liver, pancreas, testis, thymus, heart, and kidney. Increased levels are seen in hepatocellular carcinoma and pancreatic adenocarcinoma.


### 配列類似性

Contains 1 SAP domain.

### 細胞内局在

Nucleus.

## 画像



Samples: whole cell lysate (1 mg for IP, 20% of IP loaded) from HeLa cells.

Antibodies: ab115175 used for IP at 6 µg/mg lysate.

WB was then performed with a different anti-CIP29 antibody.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors