

Product datasheet

Anti-FMO1 antibody ab116251

画像数 1

製品の概要

製品名	Anti-FMO1 antibody
製品の詳細	Rabbit polyclonal to FMO1
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Rat 交差が予測される動物種: Mouse, Rabbit, Horse, Guinea pig, Cow, Cat, Dog, Human, Pig
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 40-89 (LWRFTEHVEE GRASLYNSVV SNSSKEMSCY SDFPFPEYDYP NFVPNSLFLE) of Rat FMO1 (NP_036924). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Rat muscle lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab116251** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 60 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能	This protein is involved in the oxidative metabolism of a variety of xenobiotics such as drugs and pesticides. Form I catalyzes the N-oxygenation of secondary and tertiary amines.
組織特異性	Expressed mainly in fetal liver, adult kidney and, to a lesser extent, the intestine.
配列類似性	Belongs to the FMO family.
細胞内局在	Microsome membrane. Endoplasmic reticulum membrane.

画像



Anti-FMO1 antibody (ab116251) at 1 μ g/ml + Rat muscle lysate at 10 μ g

Predicted band size: 60 kDa

12% gel

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors