


Product datasheet

Anti-KLHL12 antibody ab116252

画像数 1

製品の概要

製品名	Anti-KLHL12 antibody
製品の詳細	Rabbit polyclonal to KLHL12
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Rat 交差が予測される動物種: Mouse, Rabbit, Horse, Guinea pig, Cow, Cat, Dog, Human, Zebrafish 
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 9-58 (DIMTNTHAKS ILNSMNSLRK SNTLCDVTLR VEQKDFPAHR MLAACSDYF) of Rat C3IP1 (NP_714952.1). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Rat Brain lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituents: 98% PBS, 2% Sucrose
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab116252** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 63 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能	Serves as a substrate-specific adapter for the CUL3-based ubiquitin-protein E3 ligase complex. Negatively regulates the Wnt signaling pathway via the targeted ubiquitination and subsequent proteolysis of DVL3.
組織特異性	Highly expressed in testis and at lower levels in the submandibular salivary gland.
配列類似性	Contains 1 BTB (POZ) domain. Contains 6 Kelch repeats.
ドメイン	The BTB domain is required for interaction with CUL3.

画像



Anti-KLHL12 antibody (ab116252) at 1 µg/ml
+ Rat Brain lysate at 10 µg

Predicted band size: 63 kDa

Gel concentration: 12%.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors