


Product datasheet

Anti-Myosin IIIB antibody - N-terminal ab117247

画像数 1

製品の概要

製品名	Anti-Myosin IIIB antibody - N-terminal
製品の詳細	Rabbit polyclonal to Myosin IIIB - N-terminal
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse 
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 299-329 of Mouse Myosin IIIB (NP_796350).
ポジティブ・コントロール	NCI H460 cell line lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium azide Constituent: 99% PBS
精製度	Immunogen affinity purified
特記事項(精製)	ab117247 is purified through a protein A column, followed by peptide affinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab117247** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/100 - 1/500. Predicted molecular weight: 152 kDa.

## ターゲット情報

機能	Probable actin-based motor with a protein kinase activity.
組織特異性	Expressed in retina, kidney and testis.
配列類似性	In the N-terminal section; belongs to the protein kinase superfamily. STE Ser/Thr protein kinase family. Contains 2 IQ domains. Contains 1 myosin head-like domain. Contains 1 protein kinase domain.
細胞内局在	Cytoplasm > cytoskeleton.

## 画像



Anti-Myosin IIIB antibody - N-terminal  
(ab117247) at 1/100 dilution + NCI H460 cell  
line lysate at 35 µg

**Predicted band size:** 152 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors