

## Product datasheet

# Anti-KLHL2 antibody [8D8] ab117934

### 画像数 1

#### 製品の概要

---

製品名	Anti-KLHL2 antibody [8D8]
製品の詳細	Mouse monoclonal [8D8] to KLHL2
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment, corresponding to amino acids 1-100 and 494-593 of Human KLHL2 produced in HEK293T cells (NP_009177).
ポジティブ・コントロール	HEK293T cell lysate transfected with pCMV6-ENTRY KLHL2 cDNA.
特記事項	Stable for 12 months from date of receipt.

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 48% PBS, 1% BSA, 50% Glycerol
精製度	Protein A purified
特記事項(精製)	ab117934 is purified from Mouse ascites fluids by affinity chromatography.
ポリ/モノ	モノクローナル
クローン名	8D8
アイソタイプ	IgG2b

#### アプリケーション

---

Our [Abpromise guarantee](#) covers the use of **ab117934** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500. Predicted molecular weight: 65 kDa.

## ターゲット情報

<b>機能</b>	May play a role in organizing the actin cytoskeleton of the brain cells.
<b>組織特異性</b>	Predominantly expressed in brain.
<b>配列類似性</b>	Contains 1 BTB (POZ) domain. Contains 6 Kelch repeats.
<b>細胞内局在</b>	Cytoplasm > cytoskeleton.

## 画像



**All lanes :** Anti-KLHL2 antibody [8D8] (ab117934) at 1/500 dilution

**Lane 1 :** HEK293T cell lysate transfected with pCMV6-ENTRY control cDNA

**Lane 2 :** HEK293T cell lysate transfected with pCMV6-ENTRY KLHL2 cDNA

Lysates/proteins at 5 µg per lane.

**Predicted band size:** 65 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC205347)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors