

Product datasheet

Anti-IDH3A antibody [2E9] ab118035

画像数 2

製品の概要

製品名	Anti-IDH3A antibody [2E9]
製品の詳細	Mouse monoclonal [2E9] to IDH3A
由来種	Mouse
アプリケーション	適用あり: Flow Cyt, ICC/IF
種交差性	交差種: Human
免疫原	Full length human recombinant protein of Human IDH3A produced in HEK293T cells (NP_005521).
ポジティブ・コントロール	COS7 cells transiently transfected by pCMV6-ENTRY IDH3A, HEK293T cells transfected with pCMV6-ENTRY IDH3A

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 48% PBS, 1% BSA, 50% Glycerol
精製度	Protein G purified
特記事項(精製)	ab118035 was purified from mouse ascites fluids by affinity chromatography.
ポリ/モノ	モノクローナル
クローン名	2E9
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab118035** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt	1/100.	ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.
ICC/IF	1/100.	

ターゲット情報


配列類似性	Belongs to the isocitrate and isopropylmalate dehydrogenases family.
細胞内局在	Mitochondrion.

画像



ab118035 at 1/100 immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY IDH3A.

Immunocytochemistry/ Immunofluorescence - Anti-IDH3A antibody [2E9] (ab118035)



HEK293T cells transfected with either pCMV6-ENTRY IDH3A (Red) or empty vector control plasmid (Blue) were immunostained with 1/100 ab118035, and then analyzed by flow cytometry.

Flow Cytometry - Anti-IDH3A antibody [2E9] (ab118035)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors