

## Product datasheet

# Anti-C1orf229 antibody ab122230

### 画像数 2

#### 製品の概要

製品名	Anti-C1orf229 antibody
製品の詳細	Rabbit polyclonal to C1orf229
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB, IHC-P
種交差性	<b>交差種:</b> Human
免疫原	antigen sequence: RGAYELQGGG SQDGPGQAAV GATPTTGPGT GGEGALLGCG SGRTPTSAT WRRLLPAEV PPGAAAANFP ER, corresponding to amino acids 163-234 of Human C1orf229 (UniProt ID: Q6ZS94). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	Human colon tissue

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab122230** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/100 - 1/250.
IHC-P		1/20 - 1/50. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## ターゲット情報

## 画像




Western blot - Anti-C1orf229 antibody (ab122230)

**Lanes 2-3** : Anti-C1orf229 antibody  
(ab122230)

**Lane 1** : Marker

**Lane 2** : HEK293T cell lysate - vector only transfected

**Lane 3** : HEK293T cell lysate - overexpressing C1orf229 with a C-terminal myc-DDK tag (~3.1kDa)



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-C1orf229 antibody (ab122230)

ab122230, at a 1/30 dilution, staining C1orf229 in paraffin embedded Human colon tissue by Immunohistochemistry.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors