

Product datasheet

Anti-CDC37L1 antibody ab123006

画像数 1

製品の概要

製品名	Anti-CDC37L1 antibody
製品の詳細	Rabbit polyclonal to CDC37L1
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog
免疫原	Synthetic peptide corresponding to a region within C terminal amino acids 214-263 (EQIAHQAVVM QFIMEMAKNC NVDPRGCFRL FFQKAKAE EEE GYFEAFKNEL) of Human CDC37L1 (NP_060383). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	THP-1 cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituents: 97% PBS, 2% Sucrose
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab123006** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 39 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能	Co-chaperone that binds to numerous proteins and promotes their interaction with Hsp70 and Hsp90.
組織特異性	Expressed in brain, heart, kidney, liver, placenta and skeletal muscle.
配列類似性	Belongs to the CDC37 family.
細胞内局在	Cytoplasm.

画像



Western blot - Anti-CDC37L1 antibody (ab123006)

Anti-CDC37L1 antibody (ab123006) at 1
μg/ml + THP-1 cell lysate at 10 μg

Predicted band size: 39 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors