

## Product datasheet

# Anti-CD35 antibody [Ber-MAC DRC] ab16720

★★★★☆ 1 Abreviews

### 製品の概要

製品名	Anti-CD35 antibody [Ber-MAC DRC]
製品の詳細	Mouse monoclonal [Ber-MAC DRC] to CD35
由来種	Mouse
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Suspended human tonsillar cells.
ポジティブ・コントロール	Tonsil

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, Tissue culture supernatant
精製度	Tissue culture supernatant
ポリ/モノ	モノクローナル
クローン名	Ber-MAC DRC
アイソタイプ	IgG1

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab16720** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P	★★★★☆	

### 追加情報

IHC-P: 1/10 - 1/20.

We suggested an incubation period of 30-60 minutes at room temperature. Proteolytic treatment is required prior to immunostaining of formalin-fixed paraffin-embedded tissue sections.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

---

機能	Mediates cellular binding of particles and immune complexes that have activated complement.
組織特異性	Present on erythrocytes, leukocytes, glomerular podocytes, and splenic follicular dendritic cells.
配列類似性	Belongs to the receptors of complement activation (RCA) family. Contains 30 Sushi (CCP/SCR) domains.
細胞内局在	Membrane.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors