


## Product datasheet

# Anti-FSH beta antibody [INN-hFSH-6] ab22473

### 製品の概要

製品名	Anti-FSH beta antibody [INN-hFSH-6]
製品の詳細	Mouse monoclonal [INN-hFSH-6] to FSH beta
由来種	Mouse
特異性	ab22473 recognises the b3 epitope of FSH beta. It may be used in combination with clones INN-hFSH-117 or 132 but not hFSH beta antibodies
アプリケーション	<b>適用あり:</b> ELISA, IP, RIA
種交差性	<b>交差種:</b> Human <b>交差が予測される動物種:</b> a wide range of other species 
免疫原	Full length protein (Human): FSH 2
エピトープ	The epitope recognised by ab22473 is in the b3 region of FSH beta.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium Azide Constituents: PBS
精製度	IgG fraction
特記事項 (精製)	ab22473 is an IgG purified antibody by ion exchange chromatography from ascites.
ポリ/モノ	モノクローナル
クローン名	INN-hFSH-6
ミエローマ	x63-Ag8.653
アイソタイプ	IgG1
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab22473** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
IP		Use at an assay dependent dilution.
RIA		Use at an assay dependent dilution.

## ターゲット情報

関連性	Follicle stimulating hormone enables ovarian folliculogenesis to the antral follicle stage and is essential for Sertoli cell proliferation and maintenance of sperm quality in the testis. Members of the pituitary glycoprotein hormone family, of which FSH is one consist of a shared alpha chain and a beta chain encoded by a separate gene. The FSH beta gene encodes the beta subunit of follicle stimulating hormone.
細胞内局在	Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors