

Product datasheet

Anti-IL-5 antibody [JES1-39D10] (Biotin) ab24952

製品の概要

製品名	Anti-IL-5 antibody [JES1-39D10] (Biotin)
製品の詳細	Rat monoclonal [JES1-39D10] to IL-5 (Biotin)
由来種	Rat
標識	Biotin
アプリケーション	適用あり: ICC/IF, Flow Cyt, WB, ELISA
種交差性	交差種: Human
免疫原	The details of the immunogen for this antibody are not available.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium azide Constituent: PBS
精製度	IgG fraction
ポリ/モノ	モノクローナル
クローン名	JES1-39D10
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab24952** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. ab18445 - Rat monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 15 kDa.
ELISA		Use at an assay dependent concentration.

ターゲット情報

機能	Factor that induces terminal differentiation of late-developing B-cells to immunoglobulin secreting cells.
配列類似性	Belongs to the IL-5 family.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors