

Product datasheet

Anti-CD5 antibody [4H8E6] (Allophycocyanin) ab25280

製品の概要

製品名	Anti-CD5 antibody [4H8E6] (Allophycocyanin)
製品の詳細	Rat monoclonal [4H8E6] to CD5 (Allophycocyanin)
由来種	Rat
標識	Allophycocyanin. Ex: 645nm, Em: 660nm
アプリケーション	適用あり: WB, Flow Cyt, IP, IHC-Fr
種交差性	交差種: Mouse
免疫原	The details of the immunogen for this antibody are not available.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 16% Sucrose, PBS; Stabilizing agent
精製度	IgG fraction
ポリ/モノ	モノクローナル
クローン名	4H8E6
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab25280** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		
Flow Cyt		

アプリケーション	Abreviews	特記事項
IP		
IHC-Fr		
追加情報	<p>Flow Cyt: Use 0.2µg for 10⁶ cells.</p> <p>IHC-Fr: Use at an assay dependent dilution.</p> <p>IP: Use at an assay dependent dilution.</p> <p>WB: Use at an assay dependent dilution. Predicted molecular weight: 55 kDa.</p> <p>T cell activation. Use at an assay dependent dilution.</p> <p>Not yet tested in other applications.</p> <p>Optimal dilutions/concentrations should be determined by the end user.</p>	
ターゲット情報		
機能	May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.	
配列類似性	Contains 3 SRCR domains.	
細胞内局在	Membrane.	

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors