

Product datasheet

Anti-CD90 / Thy1 antibody [30-H12] (Phycoerythrin) ab25694

製品の概要

製品名	Anti-CD90 / Thy1 antibody [30-H12] (Phycoerythrin)
製品の詳細	Rat monoclonal [30-H12] to CD90 / Thy1 (Phycoerythrin)
由来種	Rat
標識	Phycoerythrin. Ex: 488nm, Em: 575nm
アプリケーション	適用あり: Flow Cyt, IHC-Fr
種交差性	交差種: Mouse
免疫原	Tissue, cells or virus corresponding to Mouse CD90/ Thy1. Immunogen is Mouse thymus or spleen.

特記事項

This antibody has been shown to be useful in 'in vivo' and 'in vitro' studies of T cell depletion.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.2 Preservative: 0.09% Sodium azide Constituents: 0.16% Sodium phosphate, 0.87% Sodium chloride, 0.1% Gelatin
精製度	Affinity purified
一次抗体 備考	This antibody has been shown to be useful in 'in vivo' and 'in vitro' studies of T cell depletion.
ポリ/モノ	モノクローナル
クローン名	30-H12
アイソタイプ	IgG2b
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab25694** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 0.05µg for 10 ⁶ cells.
IHC-Fr		Use at an assay dependent concentration.

ターゲット情報

機能	May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain.
配列類似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
細胞内局在	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors