

Product datasheet

Anti-V5 tag antibody [SV5-Pk1] (FITC) ab27674

2 References

製品の概要

製品名	Anti-V5 tag antibody [SV5-Pk1] (FITC)
製品の詳細	Mouse monoclonal [SV5-Pk1] to V5 tag (FITC)
由来種	Mouse
標識	FITC. Ex: 493nm, Em: 528nm
アプリケーション	適用あり: ICC/IF
免疫原	Mice were infected with the paramyxovirus SV5, Simian-Virus 5
エピトープ	Ab27674 recognises a small epitope, termed Pk, present on the P/V proteins of the paramyxovirus, SV5.
特記事項	Ab27674 has been used to detect recombinant proteins, some of which include transmembrane and secreted proteins, which have been tagged with the V5 epitope.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium Azide. Constituents: 1% BSA.
精製度	Protein G purified
特記事項 (精製)	Ab27674 was purified from tissue culture supernatant by affinity chromatography.
一次抗体 備考	Ab27674 has been used to detect recombinant proteins, some of which include transmembrane and secreted proteins, which have been tagged with the V5 epitope.
ポリ/モノ	モノクローナル
クローン名	SV5-Pk1
ミエローマ	Sp2/0-Ag14
アイソタイプ	IgG2a

アプリケーション

Our **Abpromise guarantee** covers the use of **ab27674** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		1/100.

ターゲット情報

関連性 The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKPIPPLLGLDST), although it has also been used with a shorter 9 amino acid sequence (IPNPLLGLD).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors