

Product datasheet

Anti-MSY2/YBOX2 antibody ab28643

画像数 1

製品の概要

製品名	Anti-MSY2/YBOX2 antibody
製品の詳細	Rabbit polyclonal to MSY2/YBOX2
由来種	Rabbit
アプリケーション	適用あり: ELISA, IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Cow, Dog
免疫原	A region within synthetic peptide: RRFIPRPPSV APPPMVAEIP SAGTGPGSKG ERAEDSGQRP RRWCPPPFY, corresponding to internal sequence amino acids 181-230 of Human MSY2/YBOX2 Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	HepG2 cell lysate, Human kidney tissue.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab28643** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/12500.

アプリケーション

Abreviews


特記事項

IHC-P

Use at an assay dependent concentration.

ターゲット情報

画像



Human kidney tissue stained with ab28643 at 4 µg/ml. The arrows point to the epithelial cells of the renal tubule [positive label] (400X magnification).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MSY2/YBOX2 antibody (ab28643)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors