

Product datasheet

Anti-DMRTB1 antibody ab30910

画像数 1

製品の概要

製品名	Anti-DMRTB1 antibody
製品の詳細	Rabbit polyclonal to DMRTB1
由来種	Rabbit
特異性	ab30910 recognises DMRTB1.
アプリケーション	<b>適用あり:</b> ELISA, WB
種交差性	<b>交差種:</b> Human
免疫原	A region within synthetic peptide: YGPLDLRRP MRTVGPLFT DFVRPLNINP DRALGPEYPG GSSMHPYCPF, corresponding to internal sequence amino acids 171-220 of Human DMRTB1 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	Jurkat cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab30910** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
ELISA		1/312500.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 5 µg/ml. Predicted molecular weight: 36 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## ターゲット情報

組織特異性	Testis.
配列類似性	Belongs to the DMRT family. Contains 1 DM DNA-binding domain.
細胞内局在	Nucleus.

## 画像



Anti-DMRTB1 antibody (ab30910) at 5 µg/ml  
+ Jurkat cell lysate.

### Secondary

HRP conjugated anti-Rabbit IgG  
1/50,000 - 1/100,000

**Predicted band size:** 36 kDa

**Observed band size:** 44 kDa

Suggested dilution in 5% skim milk / PBS  
buffer,

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors