

## Product datasheet

# Anti-ZNF670 antibody ab30979

### 画像数 1

#### 製品の概要

製品名	Anti-ZNF670 antibody
製品の詳細	Rabbit polyclonal to ZNF670
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human
免疫原	A region within synthetic peptide: EIKEGSQYGE TFSQDSNLNL NKKVSTGVKP CECSVCGKVF ICHSALHRHI, corresponding to N terminal amino acids 73-122 of Human ZNF670 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	HepG2 cell lysate

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab30979** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1.25 µg/ml. Detects a band of approximately 40 kDa (predicted molecular weight: 45 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

アプリ ケーション	Abreviews	特記事項
ELISA		1/312500.

## ターゲット情報

関連性	The function of zinc finger protein 670 is currently unknown. The gene is located on chromosome 1.
細胞内局在	Nuclear

## 画像



Anti-ZNF670 antibody (ab30979) at 1.25 µg/ml (in 5% skim milk / PBS buffer) + HepG2 cell lysate

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size:** 45 kDa

**Observed band size:** 40 kDa

[why is the actual band size different from the predicted?](#)

gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors