

Product datasheet

Anti-ZNF 559 antibody ab36187

画像数 1

製品の概要

製品名	Anti-ZNF 559 antibody
製品の詳細	Rabbit polyclonal to ZNF 559
由来種	Rabbit
アプリケーション	適用あり: ELISA, WB
種交差性	交差種: Human
免疫原	Synthetic peptide: QEKPYKCSDC EKGLPSSSHL RECVRIYGGE RPYTHKEYVE TFSHSTALFV, corresponding to a region within N terminal amino acids 181-230 of Human ZNF 559 Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Jurkat cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS, pH 7.4
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab36187** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項

アプリ ケーション	Abreviews	特記事項
ELISA		1/1.5625e+006. ELISA titre using peptide based assay.
WB		Use a concentration of 5 µg/ml. Detects a band of approximately 59 kDa (predicted molecular weight: 62 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能	May be involved in transcriptional regulation.
配列類似性	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 13 C2H2-type zinc fingers. Contains 1 KRAB domain.
細胞内局在	Nucleus.

画像



Anti-ZNF 559 antibody (ab36187) at 2.5 µg/ml
+ Jurkat cell lysate at 10 µg

Secondary

HRP conjugated anti-rabbit IgG (diluted to
1/50,000-1/100,000)

Predicted band size: 62 kDa

Observed band size: 59 kDa

Suggested dilution of ab36187 at 5.0µg/ml in
5% skim milk / PBS buffer

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors