

Product datasheet

Anti-PKC delta antibody ab4142

1 References

製品の概要

製品名	Anti-PKC delta antibody
製品の詳細	Rabbit polyclonal to PKC delta
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA, IHC-Fr, ICC
種交差性	交差種: Mouse, Rat 交差が予測される動物種: Mammals
免疫原	Synthetic peptide: SFVNPKYEQFLE , corresponding to amino acids 662-673 of Rat PKC delta Run BLAST with Run BLAST with
ポジティブ・コントロール	rat cerebellar lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.05% Sodium Azide
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab4142** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		

アプリケーション	Abreviews	特記事項
ELISA		
IHC-Fr		
ICC		
追加情報	<p>ELISA: Use at an assay dependent dilution. ICC: 1/250. IHC-Fr: Use at an assay dependent dilution. WB: 1/250. Predicted molecular weight: 81 kDa. Use lysis buffer made of 20 mM Tris-HCl, pH 7.5, 0.25 M sucrose, 10 mM EGTA, 1 mM EDTA, 10 mM β-mercaptoethanol, 1 mM PMSF, 0.02% leupeptin, and 0.1% Triton X-100. Incubate the homogenate for 1 h at 4°C and centrifuge at 100,000 g for 30 min at 4°C. (see Wetsel et al reference).</p> <p>Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.</p>	
ターゲット情報		
機能	<p>This is calcium-independent, phospholipid-dependent, serine- and threonine-specific enzyme. PKC is activated by diacylglycerol which in turn phosphorylates a range of cellular proteins. PKC also serves as the receptor for phorbol esters, a class of tumor promoters. May play a role in antigen-dependent control of B-cell function. Phosphorylates MUC1 in the C-terminal and regulates the interaction between MUC1 and beta-catenin.</p>	
配列類似性	<p>Belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family. PKC subfamily. Contains 1 AGC-kinase C-terminal domain. Contains 1 C2 domain. Contains 2 phorbol-ester/DAG-type zinc fingers. Contains 1 protein kinase domain.</p>	
ドメイン	<p>The C1 domain, containing the phorbol ester/DAG-type region 1 (C1A) and 2 (C1B), is the diacylglycerol sensor. The C2 domain is a non-calcium binding domain. It binds proteins containing phosphotyrosine in a sequence-specific manner.</p>	
翻訳後修飾	<p>Phosphorylated on Thr-507, within the activation loop. Autophosphorylated and/or phosphorylated. Although the Thr-507 phosphorylation occurs it is not a prerequisite for enzymatic activity.</p>	
細胞内局在	<p>Cytoplasm. Membrane.</p>	

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors