

Product datasheet

Anti-IL17 antibody [TC11-18H10] ab47794

1 References

製品の概要

製品名	Anti-IL17 antibody [TC11-18H10]
製品の詳細	Rat monoclonal [TC11-18H10] to IL17
由来種	Rat
アプリケーション	適用あり: ELISA
種交差性	交差種: Mouse
免疫原	The details of the immunogen for this antibody are not available.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	モノクローナル
クローン名	TC11-18H10
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab47794** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.

ターゲット情報

機能 Ligand for IL17RA and IL17RC (PubMed:17911633). The heterodimer formed by IL17A and

IL17F is a ligand for the heterodimeric complex formed by IL17RA and IL17RC (PubMed:18684971). Involved in inducing stromal cells to produce proinflammatory and hematopoietic cytokines (PubMed:8676080).

組織特異性

Restricted to activated memory T-cells.

配列類似性

Belongs to the IL-17 family.

翻訳後修飾

Found both in glycosylated and nonglycosylated forms.

細胞内局在

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors