

Product datasheet

Anti-CPXCR1 antibody [2205C1] ab51346

画像数 1

製品の概要

製品名	Anti-CPXCR1 antibody [2205C1]
製品の詳細	Mouse monoclonal [2205C1] to CPXCR1
由来種	Mouse
アプリケーション	適用あり: WB, Dot blot
種交差性	交差種: Human
免疫原	Recombinant human CPXCR1

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS (8mM Sodium phosphate, 3mM Potassium chloride, 140mM Sodium chloride, 1.5mM Potassium phosphate, pH 7.4).
精製度	Protein G purified
特記事項 (精製)	Ab51346 was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG depleted (approximately 95%) fetal bovine serum.
ポリ/モノ	モノクローナル
クローン名	2205C1
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab51346** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution. Predicted molecular weight: 35 kDa.

アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent dilution.

ターゲット情報

組織特異性 Expressed in a variety of fetal tissues.

画像



Anti-CPXCR1 antibody [2205C1] (ab51346) + Recombinant CPXCR1

Predicted band size: 35 kDa

Observed band size: 44 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors