

Product datasheet

Anti-Salivary alpha amylase antibody ab54765

1 References [画像数 3](#)

製品の概要

製品名	Anti-Salivary alpha amylase antibody
製品の詳細	Mouse monoclonal to Salivary alpha amylase
由来種	Mouse
アプリケーション	適用あり: IHC-P, WB
種交差性	交差種: Human, Recombinant fragment
免疫原	Recombinant fragment corresponding to Human Salivary alpha amylase aa 172-246. Sequence: VRDCRLSGLLDLALGKDYVRSKIAEYMNHLIDIGVAGFRIDASKHMWPGD IKAILDKLHNLNSNWFPEGSKPFI

 [Run BLAST with](#)

 [Run BLAST with](#)

特記事項

Abcam is committed to meeting high standards of ethical manufacturing and has decided to discontinue this product by June 2019 as it has been generated by the ascites method. We are sorry for any inconvenience this may cause. We would recommend antibody [ab201450](#) as a replacement.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab54765** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 3 µg/ml.
WB		Use a concentration of 1 - 5 µg/ml. This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

ターゲット情報


配列類似性	Belongs to the glycosyl hydrolase 13 family.
細胞内局在	Secreted.

画像




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Salivary alpha amylase antibody (ab54765)

alpha Amylase antibody (ab54765) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human salivary gland.



Western blot - Anti-Salivary alpha amylase antibody (ab54765)

Western blot against tagged recombinant protein immunogen using ab54765 alpha Amylase antibody at 1ug/ml. Predicted band size of immunogen is 34 kDa



Anti-Salivary alpha amylase antibody
(ab54765) + Human pancreas

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors