

## Product datasheet

# Anti-SFRS14 antibody ab56575

画像数 2

### 製品の概要

製品名	Anti-SFRS14 antibody
製品の詳細	Mouse monoclonal to SFRS14
由来種	Mouse
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Recombinant fragment: PQRADHRVVG TIDQLVKRVI EGSLSPKERT LLKEDPAYWF LSDENSLEYK YYKLKLAEMQ RMsENLRGAD QKPTSADCAV RAMLYSRAVR NLKKKL, corresponding to amino acids 579-675 of Human SFRS14 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab56575** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 120 kDa.

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 3 µg/ml.

## ターゲット情報


機能	May play a role in mRNA splicing.
組織特異性	Detected in adult testis, and in fetal brain and kidney.
配列類似性	Contains 1 G-patch domain. Contains 2 SURP motif repeats.
細胞内局在	Nucleus.

## 画像



SFRS14 antibody (ab56575) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human testis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SFRS14 antibody (ab56575)



SFRS14 antibody (ab56575) at 1ug/lane + HeLa cell lysate at 25ug/lane.

Western blot - Anti-SFRS14 antibody (ab56575)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors