

Product datasheet

Anti-CTNNBL1 antibody ab57617

画像数 2

製品の概要

製品名	Anti-CTNNBL1 antibody
製品の詳細	Mouse monoclonal to CTNNBL1
由来種	Mouse
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Recombinant fragment, corresponding to amino acids 210-311 of Human CTNNBL1

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG3
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab57617** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 65 kDa.
IHC-P		Use a concentration of 3 µg/ml.

ターゲット情報


機能	Induces apoptosis in CHO cells.
組織特異性	Widely expressed with highest levels in skeletal muscle, placenta, heart, spleen, testis and thyroid.
細胞内局在	Nucleus.

画像



CTNNBL1 antibody (ab57617) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human testis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CTNNBL1 antibody (ab57617)



CTNNBL1 antibody (ab57617) at 1ug/lane + HeLa cell lysate at 25ug/lane.

Western blot - Anti-CTNNBL1 antibody (ab57617)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors