

## Product datasheet

# Anti-EFHD1 antibody ab57818

画像数 3

### 製品の概要

製品名	Anti-EFHD1 antibody
製品の詳細	Mouse monoclonal to EFHD1
由来種	Mouse
アプリケーション	適用あり: WB, ICC/IF
種交差性	交差種: Human
免疫原	Recombinant full length protein: GLMALAKLSE IDVALEGVKG AKNFFEAKVQ ALSSASKFEA ELKAEQDERK REEEERRLRQ AAFQKLFK N, corresponding to amino acids 168-238 of Human EFHD1 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリモノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab57818** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml.

## アプリケーション

## Abreviews

## 特記事項

ICC/IF


Use a concentration of 10 µg/ml.

## ターゲット情報

### 配列類似性


Contains 2 EF-hand domains.

## 画像




Immunohistochemical analysis of Formalin/PFA-fixed paraffin-embedded human kidney sections labeling EFHD1 using ab57818 at antibody concentration 3 µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-EFHD1 antibody (ab57818)



Immunofluorescence analysis of HeLa cell lysates labeling EFHD1 with ab57818 at 10 µg/ml.

Immunocytochemistry/ Immunofluorescence - Anti-EFHD1 antibody (ab57818)



Western blot against tagged recombinant protein immunogen using ab57818 EFHD1 antibody at 1ug/ml. Predicted band size of immunogen is 34 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors