

Product datasheet

Anti-MEX3B antibody ab57910

画像数 1

製品の概要

製品名	Anti-MEX3B antibody
製品の詳細	Mouse monoclonal to MEX3B
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment: RKKSVNMTEC VPVPSSSEHVA EIVGRQGCKI KALRAKTNTY IKTPVRGEEP VFVVTGRKED VAMARREIIS AAHFMSMIRA SRNKNTALNG AVPGPPNLPG , corresponding to amino acids 60-160 of Human RKHD3 Run BLAST with ExPASy Run BLAST with NCBI

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG2b
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab57910** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項


アプリケーション	Abreviews	特記事項
----------	-----------	------

WB Use a concentration of 1 - 5 µg/ml.
This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

ターゲット情報

機能	RNA-binding protein. May be involved in post-transcriptional regulatory mechanisms.
組織特異性	Highest levels found in fetal brain and testis. Detected in the adult intestinal epithelium, specifically in goblet cell at protein level.
配列類似性	Contains 2 KH domains. Contains 1 RING-type zinc finger.
ドメイン	Binds RNA through its KH domains.
翻訳後修飾	Phosphorylated.
細胞内局在	Nucleus. Cytoplasm. Cytoplasm > P-body. Predominantly expressed in the cytoplasm and shuttles between the cytoplasm and the nucleus through the CRM1 export pathway.

画像



Western blot against tagged recombinant protein immunogen using ab57910 RKHD3 antibody at 1ug/ml. Predicted band size of immunogen is 37 kDa

Western blot - Anti-RKHD3 antibody (ab57910)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors