

Product datasheet

Anti-SPRR3 antibody ab58233

2 References [画像数 2](#)

製品の概要

製品名	Anti-SPRR3 antibody
製品の詳細	Mouse monoclonal to SPRR3
由来種	Mouse
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human
免疫原	Recombinant full length protein, corresponding to amino acids 1-162 of Human SPRR3

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None PBS, pH 7.2
精製度	Protein G purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG3
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab58233** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml.
IHC-P		Use a concentration of 3 µg/ml.

ターゲット情報


機能	Cross-linked envelope protein of keratinocytes.
配列類似性	Belongs to the cornifin (SPRR) family.
細胞内局在	Cytoplasm.

画像



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SPRR3 antibody (ab58233)

SPRR3 antibody (ab58233) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human tonsil.



Western blot - Anti-SPRR3 antibody (ab58233)

Western blot against tagged recombinant protein immunogen using ab58233 SPRR3 antibody at 1ug/ml. Predicted band size of immunogen is 44 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors